



Reliability, accuracy, and precision are essential for any quality pH measurement. In order to guarantee the best quality measurements, sensors must be suited for the sample, easy to use and designed to perform.

The new InLab® sensor portfolio combines high-quality materials with intuitive design and offer a blend of new and tried-and-trusted technologies, including:

- ISM® technology for automatic sensor recognition
- Innovative new wetting cap for better storage and easy handling
- Optimised sensor build for demanding and challenging samples

Learn more how the new InLab sensor portfolio can help you improve the quality of your pH measurements by downloading the 'Origin of Precision' sensor poster - simply visit [www.mt.com/uk-ph-sensor](http://www.mt.com/uk-ph-sensor)

For more information call 0116 234 5095 or alternatively email [enquire.mtuk@mt.com](mailto:enquire.mtuk@mt.com) .